

January  
1955

U. S. DEPARTMENT OF COMMERCE  
NATIONAL BUREAU OF STANDARDS  
WASHINGTON 25, D. C.

Letter  
Circular  
LC-1017

STANDARDS FOR CHECKING THE  
CALIBRATION OF SPECTROPHOTOMETERS  
(200 to 1000  $m\mu$ )

Contents

1. Introduction.
  2. Checking the wavelength scale.
    - 2.1. Non-recording spectrophotometers.
    - 2.2. General Electric recording spectrophotometers.
  3. Checking the photometric scale.
    - 3.1. Glass standards of spectral transmittance.
    - 3.2. Solution standards of spectral transmittancy.
  4. Reference standards of spectral directional reflectance calibrated relative to freshly prepared magnesium oxide.
    - 4.1. For the General Electric recording spectrophotometer.
    - 4.2. For the Beckman quartz (model DU, non-recording) spectrophotometer.
  5. References.
1. Introduction.

In continuation of a type of activity carried on for many years at the National Bureau of Standards, there is described in this letter circular the various types of standards that are available for issuance by the Bureau for the purpose of checking or maintaining the over-all reliability of spectrophotometers in the ultraviolet, visible and near-infrared regions of the spectrum. Certain other information of similar purpose is also included.

Experience has shown that spectrophotometers can easily get out of adjustment. Although repeated trials may give the same values over and over again, indicating high sensitivity and precision, and the instrument may otherwise appear to be functioning perfectly, gross errors in wavelength may nevertheless be rendering the values obtained highly unreliable. Other causes of error may likewise be present and unsuspected, particularly with the photoelectric spectrophotometers now comprising so large a percentage of the total in use. The use of the various standards described herein has been found of considerable assistance in detecting and eliminating errors that would otherwise be present, or in confirming that the instrument is in fact giving reliable results.

Reference is made in this letter circular to the NBS test fee schedules currently in effect. Those relating to spectrophotometry are designated as 202.185, Spectrophotometric Standards, and 202.186, Spectrophotometric Measurements. These schedules are obtainable from the National Bureau of Standards on demand.

## 2. Checking the Wavelength Scale.

Most of the present-day spectrophotometers have a direct-reading wavelength scale; that is, the scale, instead of being divided in uniform linear or circular measure, is divided and engraved directly in millimicrons. This greatly facilitates setting the instrument at any desired wavelength. The accuracy of many of these direct-reading wavelength scales is remarkably good, when put in the best average adjustment, considering the difficulties of quantity production of such scales. When so adjusted, it is not uncommon to find them in error by not more than 1  $m\mu$  throughout the ultraviolet and visible spectrum. However, if one wishes the uncertainties in his wavelength settings to be of the order of 0.1  $m\mu$ , a careful check of these direct-reading scales is necessary.

Certain sources and wavelengths that have proved especially suitable for the calibration of various types of spectrophotometers are listed in table 1, which is similar in scope and purpose to table 1 of Circular 484 (1)\*.

---

\*Numbers in parentheses refer to the References, Section 5.

---

All values of table 1 are consistent with those published in the M.I.T. wavelength tables (2). Which of these sources to use, and which wavelengths of the several sources are the most suitable, will depend on the type of instrument. Furthermore, the procedure will vary importantly depending on whether the instrument is a non-recording or a recording spectrophotometer.

Table 1. Sources and Wavelengths Suitable for the Calibration of Spectrophotometers.

Mercury Arc in Quartz  
(Same Wavelengths in  
glass above 300 mμ)

Aluminum  
Spark in Air

Helium Discharge  
Tube in Glass

<u>Wavelength</u> mμ	<u>Footnote</u>	<u>Wavelength</u>	<u>Wavelength</u>	
205.29		216.88	318.77	
222.47		217.40	361.36	
223.41			363.42	
225.88		220.46	370.50	
230.21		221.00	381.96	
232.32			388.86*	
235.25		226.35	396.47	
237.83*		226.91	402.62	
239.94} *	1		412.08	
239.97} *		236.71	414.38	
244.69		237.21	438.79	
246.41		237.31	443.75	
248.20		237.34	447.15*	
248.27} *	2	237.84	471.31*	
248.38			492.19	
253.48} *	3	256.80	501.57*	
253.65} *		257.51	504.77	
257.63		257.54	587.56*	
260.32			667.81*	
265.20		263.16	706.52*	
265.37} *	4		728.13	
265.51} *		265.25	1083.0 *	
269.95		266.04		
275.28*				
275.97		281.62		Neon Discharge Tube
280.35} *	5			
280.45} *		308.22	Wave-	Relative
284.78		309.27	length	Intensity
289.36*				
292.54		358.69	585.25	5
296.73*			588.19	4
302.15		394.40	594.48	8
302.35} *	6	396.15	597.55	2
302.56} *			603.00	2
302.75} *			607.43	8
312.57*			609.62	13
313.15} *	7		614.31	25
313.18} *			616.36	6
334.15			621.73	4
349.28				

\* These lines have been found most useful on the Beckman DU spectrophotometer.

Table 1. (Continued)

Mercury Arc in Quartz (Same wavelengths in glass above 300 mμ.)		Hydrogen Sodium Cesium Arcs	Neon Discharge Tube
Wavelength mμ	Footnote		Relative Intensity
365.01	8	H 434.05	11
365.48*		H 486.13	4
366.29		Na 589.00	20
366.33		Na 589.59	23
390.64		H 656.28	100
398.40	9	Cs 852.11	39
400.63		Cs 894.35	8
404.66*			12
407.78*			23
435.83*			14
491.60			23
546.07*			2
576.96			45
579.07*			--
623.44			5
671.62			17
690.72			4
1014.0*			--
1128.7			--
			--

1. A value of 239.95 is recommended for the unresolved pair.
2. A value of 248.3 is recommended when the 3 lines are unresolved.
3. The intensity of 253.48 is negligible compared to that of 253.65. The latter value should be used when the lines are unresolved.
4. The 265.20 line is somewhat stronger than the others and a value of 265.3 is recommended when the three lines are unresolved.
5. These two lines are of approximately the same intensity and a value of 280.40 is recommended for the unresolved pair.
6. The two shorter lines are considerably stronger than the other two. It is probable that a value of 302.25 should be used for the unresolved lines.
7. A value of 313.16 is recommended for the unresolved pair.
8. With the arc used on the Beckman DU spectrophotometer the ratio of intensities for 365.01 : 365.48 : 366.33 is 100 : 48 : 36, approximately. The intensity of the 366.29 line appears negligible relative to that of 366.33.
9. These two lines are of approximately the same intensity and a value of 578.0 is recommended for the unresolved pair.

## 2.1 Non-Recording Spectrophotometers.

The best procedure for checking the wavelength scale of a non-recording spectrophotometer is by direct use of a source of radiant energy having spectral lines of suitable intensity and adequately spaced throughout the spectral range of interest. Various sources are available and can be recommended for such purpose. How many sources, or how many wavelengths, to use in such a calibration depends, of course, on the desires of the individual investigator.

In this connection it should be noted that the number of significant figures of importance in spectrophotometry (including "absorption spectroscopy") is of a different order of magnitude than that used in emission spectroscopy or in standard wavelength tables. In the visible spectrum with the usual type of spectrophotometer it seems impossible to maintain the wavelength calibration with uncertainties less than about 0.1  $m\mu$ . While the uncertainty may be less in the ultraviolet with a prism instrument, there seems no purpose served in giving standard wavelengths to better than 0.01  $m\mu$  for spectrophotometric calibration.

Two suitable sources for wavelength calibration are the mercury lamp and the helium lamp. A mercury lamp in a quartz envelope is by far the best single source for wavelength calibration from 205 to 1014  $m\mu$ . A mercury lamp in a glass envelope provides the same spectral lines except that below about 300  $m\mu$  they are not transmitted by the glass envelope.

The helium lines are especially well placed for wavelength calibration in the visible spectrum, and the strong lines at 388 and 1083  $m\mu$  are also often very useful. Many other sources, flame or arc, are available for visual wavelength calibration (2, 3) but most of these are too unstable for accurate calibration with a photoelectric detector.

These same sources and many others are also useful for the wavelength calibration of spectrographs used in photographic spectrophotometry. Between 200 and 400  $m\mu$  the series of doublets obtained from the aluminum spark in air is very useful because they are so readily recognized.

Not all of the lines for any of the sources are given in table 1 but only those that are considered especially suitable for the purpose. Furthermore, not even all of those listed for any one source may be suitable for any one particular instrument. The mercury arc in quartz is an example. All of the lines listed (and still others) can be used for wavelength

calibration of a photographic spectrophotometer over the range of sensitivity of the plate used. The lines from 404.7 to 690.7  $m\mu$  can be used for visual calibration of a spectrophotometer. But not all of the lines are suitable for calibration of a photoelectric instrument, and those that prove adequate will depend on the sensitivity and slit widths characteristic of any particular instrument. One must be very careful that other lines are not included, in addition to the one on which the settings are supposedly being made, of sufficient intensity to affect the wavelength setting.

Special attention should perhaps be called to the use of a cesium arc at 852.1 and 894.3  $m\mu$  (4). From table 1 it is apparent that there are few suitable lines between 706.5 and 1014.0  $m\mu$ , particularly from steady sources necessary or desirable in the calibration of photoelectric spectrophotometers. The neon discharge tube gives many lines between 750 and 1000  $m\mu$  (2) but these have not been found satisfactory in the calibration of photoelectric spectrophotometers. In the orange and red the neon lines are useful for visual calibration and many of these can be used to calibrate photoelectric spectrophotometers (5) if the sensitivity is such that very narrow slits can be used. The relative intensities (6) given in table 1 will help in case of overlapping.

The best technique to use in wavelength calibration of non-recording spectrophotometers, given a suitable source, will vary from instrument to instrument and method to method. A few general principles can be given here, however.

In photographic spectrophotometry it usually is sufficient to photograph a known spectrum at the top and bottom of the plate, unless the source used for the absorption spectra itself carries such known reference lines. A few of these reference lines will then serve to correlate that particular plate with whatever complete calibration curve has previously been established by more extensive measurements with the various sources.

On visual and photoelectric non-recording spectrophotometers, it usually is necessary, for highest precision, to have a basic reference line to which all of the other wavelengths are compared by direct check. At the Bureau the Hg yellow lines have proved most suitable for the König-Martens visual spectrophotometer (7). At the slit widths used the overlapping of the two lines gives a central brighter "line" taken as 578.0  $m\mu$  with a luminous background against which the slit jaws are readily seen. A luminous background, or slight illumination of the ocular slit, always facilitates calibration when an eyepiece is used. Visual calibration without an eyepiece is usually less precise unless very narrow slits are used.

Two techniques have been used at the Bureau in the calibration of non-recording photoelectric spectrophotometers. On the Gibson spectrophotometer (8) the slits are always 0.1 mm wide or greater and the most reliable calibration is obtained by plotting galvanometer deflections at closely adjacent wavelengths. The most probable value for the wavelength reading is given by the intersection of the two straight lines resulting from a plot of the data for any given line, the correction being given by the difference between this value and the true wavelength. This is illustrated in reference (1).

On the Beckman DU spectrophotometer the same method has been used (9), but at the Bureau it has seemed preferable and is much more rapid, to calibrate with a narrow slit and record the wavelength dial reading for the maximum left deflection of the galvanometer needle as the wavelength dial is slowly turned. The most suitable reference line on two of the Bureau's instruments has proved to be the Hg green line at 546.07  $m\mu$  (5).

## 2.2. General Electric Recording Spectrophotometers.

The initial wavelength calibration of a recording spectrophotometer, such as the manufacturer must carry out in connection with cutting his cams or preparing his reading scale, is not here considered, but only the check of such a calibration by the user of the instrument.

Such a user can, of course, follow the procedure prescribed above for checking the wavelength calibration of non-recording spectrophotometers. However, there are two important reasons for following a different procedure for recording spectrophotometers. For such instruments it is desirable to have a calibration that is made with the instrument operating. It is further desirable in most kinds of work to have this calibration appear on the graph sheet so that difficulties connected with positioning of the sheet, expansion or contraction of the paper with humidity or temperature, or instrumental variations can be eliminated.

Wavelength calibrations of this kind can be made if a material is available having a number of strong and narrow absorption or transmission bands suitably spaced over the spectral range of interest. Two materials have been used or suggested for this purpose: (a) Didymium glasses have been used for many years at the National Bureau of Standards (10), (b) quartz-Polaroid combinations have been proposed (11) and may prove useful for such work.

The use of a didymium glass in this manner would not in general be accurate unless it is calibrated at nearly the same slit widths as are to be used. Most of the absorption bands that are usable for the purpose are multiple bands and the wavelengths of maximum absorption often depend on the slit widths. This has been illustrated in previous publications (1, 10).

The NBS didymium glass standards were carefully calibrated by point-by-point measurements on the König-Martens visual and Gibson photoelectric spectrophotometers with slit widths approximating the 4, 8, 10 and 20 millimicron slits used on the NBS General Electric spectrophotometers. Some of these values have been published (12).

The most suitable didymium glass for the purpose, considering type of curve and availability, is a Corning 5120 glass of 3.0 mm thickness. While it is not known how much the wavelengths of maximum absorption of this 5120 glass might vary from melt to melt, glasses from at least three melts have been measured and there has never been any certain variation among the samples tested. For much work it is probably safe to use the values given in table 2.

Table 2. Wavelengths of maximum absorption for Corning 5120 glasses of 3.0 mm thickness as obtained at the National Bureau of Standards for the slit widths indicated.

Wavelength of Maximum Absorption (m $\mu$ )	Approximate Spectrum Interval Transmitted by Slits (m $\mu$ )
441.0	10
475.5	10
528.7	10
585.0	10
684.8	10
743.5	10
745.5	20
808.	20
883.	20
1067.	20

For those who wish greater certainty, however, the Bureau has obtained a supply of Corning 5120 glass in 2-inch polished squares and of 3.0 mm thickness. These are measured and the values reported in accordance with NBS test fee schedule 202.185, items d to f. The measurements consist of recording a curve of the test glass on the same sheet as the curve of the NBS standard glass and deriving values of the wavelengths of minimum transmittance of the former relative to those of the latter. The over-all uncertainties of the values so reported are considered to be not greater than  $\pm 1 \text{ m}\mu$  from 441.0  $\text{m}\mu$  to 743.5  $\text{m}\mu$ , and not greater than  $\pm 2 \text{ m}\mu$  from 745 to 1067  $\text{m}\mu$ .

Methods of use of a calibrated didymium glass on a G. E. recording spectrophotometer are described in references (1), (10), and (12).

While the use of a didymium glass for checking the wavelength calibration of G. E. recording spectrophotometers is highly recommended, as noted above, there are two other uses of the didymium glass which are not recommended. First, the didymium glass is not well suited for checking the photometric scale of any spectrophotometer, recording or non-recording. Transmittances at the peaks of the absorption bands are too dependent on slit widths, and transmittances on the steep parts of the curve are too dependent on slight wavelength errors, both as illustrated in Fig. 1 of reference (10) and in Fig. 8 of reference (1). Second, the use of a didymium glass to check the wavelength calibration of a non-recording spectrophotometer is considered much inferior from the standpoints of time, convenience and reliability to the direct use of line sources as described in Sec. 2.1. The National Bureau of Standards has consistently refused to accept didymium glasses for calibration for either of these two purposes.

### 3. Checking the Photometric Scale.

A check of the photometric scale of a spectrophotometer independent of all other sources of error is difficult or impossible to make on most spectrophotometers. Useful for this purpose would be a series of samples whose respective transmittances do not vary with wavelength, which will not displace the beam, which do not reflect strongly, and for which the transmittances can be independently determined with high accuracy. However, no such glasses or other materials are available.

On a visual spectrophotometer such a check is possible by means of rapidly rotating sectors. If such sectors are properly made, the angle of the total opening relative to  $360^\circ$ , that is, the effective transmittance of the sector, can be measured on a circular dividing engine with uncertainties only in the fifth decimal place. Of course, the use of such rotating sectors is based on the validity of Talbot's law. Throughout the spectrum, no certain deviations from this relation are known. As a matter of fact, average values obtained over a period of years in measuring the "transmittances" of accurately calibrated rotating sectors on the NBS König-Martens spectrophotometer prove both the validity of Talbot's law at various wavelengths and the reliability of the instrument over most of the photometric scale, or else there is a remarkable balancing of errors.

### 3.1 Glass Standards of Spectral Transmittance.

Usually on photoelectric spectrophotometers the use of rotating sectors to check the reliability of the photometric scale either is impossible or is attended with too much uncertainty for one reason or another.

Accordingly, shortly after the advent of commercial photoelectric spectrophotometers, the National Bureau of Standards instituted the service of issuing glass standards of spectral transmittance (13). To date, about 350 of these filters have been issued with accompanying certificates. Further information about these filters, particularly with reference to their permanence, may be found in a recent publication (14).

All of the filters offered for this purpose have transmittances varying through the visible spectrum from high (0.75 or more) to low (0.06 or less). In one sense these might be considered inferior to strictly neutral filters in that a deviation from the true value may be attributed to other causes than inaccuracy of the photometric scale. On the other hand, they are superior to the neutral filters in detecting stray-energy, slit-width, and gross wavelength errors.

The four types of glass filters are designated as "carbon yellow", "cobalt blue", "copper green", and "selenium orange". They are approximately 30 mm in diameter and 2, 3, 2, and 2 mm thick, respectively. The copper green glasses are also obtainable as 25 mm squares, 2 mm thick. They are issued in accordance with NBS test fee schedule 202.185, item a. The transmittances reported are for  $25^\circ\text{C}$ .

and are usually obtained by measuring wavelength by wavelength on the Beckman DU spectrophotometer, the ratio of transmittance of test glass to standard. The standards have been calibrated (5) by extensive measurements on the Beckman DU, König-Martens and Gibson spectrophotometers, from which also are derived the effects of temperature change on the transmittances and the major part of the uncertainties reported for the values.

To give the reader a better idea of the transmittances he may expect on the standards issued, there are given in table 3 the transmittances of the respective NBS standards at the wavelengths used and reported. The transmittances of the standards issued will not in general be identical with those of table 3 but will not be greatly different from them.

Transmittances of these filters at wavelengths other than those given will be determined on request from 365  $m\mu$  to 1000  $m\mu$  in accordance with item c of NBS test fee schedule 202.185. Values will be obtained for temperatures of 25°C. The effect of change of temperature has not been determined at the Bureau for these glasses outside the range from 390 to 750  $m\mu$ . It is known, however, that for all four types of glass the temperature effects are very small from 750 to 1000  $m\mu$ , and are probably negligible for the usual room temperature variations. On the other hand, temperature effects are always large for these kinds of glass when the transmittance curve is decreasing rapidly towards shorter wavelengths (1, Figs. 29 and 30), so that increasingly large temperature effects may be expected for these filters in the ultraviolet.

### 3.2. Solution Standards of Spectral Transmittancy.

The photometric scale of spectrophotometers may be checked by means of solutions of known spectral transmittancies, instead of by means of the standard glasses; if one prefers. For this purpose the following are recommended:

1. The aqueous solutions of copper sulphate and cobalt ammonium sulphate used in the series of filters developed at the National Bureau of Standards for reproducing the colors of sunlight and daylight and for the determination of color temperatures (15).

Table 3. Values of Spectral Transmittance of NBS Glass Standards for Checking the Photometric Scale of Spectrophotometers. The transmittances of glasses issued by the Bureau will not, in general, be identical with those of table 3 but will not be greatly different from them.

Wavelength mμ	Transmittance for glasses designated as			
	carbon yellow	cobalt blue	copper green	selenium orange
390	0.025	0.895	0.862	
404.7	.020	.884	.877	
420	.019	----	----	
435.8	.0240	.806	.893	
471.3	.081	.612	.894	
491.6	----	.344	----	
501.6	.208	.245	.859	
520	.316	.091	----	
530	.379	----	.760	
540	----	.0308	----	
546.1	.479	.0335	.671	
560	.557	.064	.585	
570	----	----	----	
578	.636	.0273	.473	0.0042
587.6	.668	----	----	.118
600	.699	.0074	.350	.55
620	.731	.0100	.256	.852
640	.747	.0074	.187	.904
660	.754	----	----	.914
667.8	----	.034	.122	----
680	----	.14	----	----
690	.755	.34	----	.919
706.5	-----*	.713	----	----
710	----	----	.074	----
720	.748	.845	----	.918
750	.730	.901	.057	.917

Table 4. Values of Spectral Absorbancy,  $A_s$ , and Transmittancy,  $T_s$  for Standard Copper Sulphate Solution as Specified, Unfiltered, Thickness 10.00 mm. Temperature 25°C, Having the Following Composition:

Copper Sulphate ( $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ )\* 20.000 grams  
 Sulphuric Acid (specific gravity 1.835) 10.0 cc  
 Water (distilled) to make 1000. cc

Wavelength $m\mu$	$A_s$	$T_s^{**}$	Wavelength $m\mu$	$A_s$	$T_s^{**}$
350	0.0090	0.979	600	0.0680	0.855
60	.0063	.986	10	.0885	.816
70	.0046	.989	20	.1125	.772
80	.0035	.992	30	.143	.719
90	.0028	.994	40	.180	.661
400	.0023	.995	650	.224	.597
10	.0019	.996	60	.274	.532
20	.0016	.996	70	.332	.466
30	.0014	.997	80	.392	.406
40	.0012	.997	90	.459	.348
450	.0011	.997	700	.527	.297
60	.0011	.997	10	.592	.256
70	.0012	.997	20	.656	.221
80	.0014	.997	30	.715	.193
90	.0018	.996	40	.768	.171
500	.0026	.994	750	.817	.152
10	.0038	.991	Hg 404.7	.0021	.995
20	.0055	.987	Hg 435.8	.0013	.997
30	.0079	.982	Hg 491.6	.0019	.996
40	.0111	.975	Hg 501.6	.0028	.994
550	.0155	.965	Hg 546.1	.0135	.969
60	.0216	.951	Hg 578.0	.0368	.919
70	.0292	.935	He 587.6	.0487	.894
80	.0390	.914	He 667.8	.319	.480
90	.0518	.888			

\*Analysis showed the copper sulphate to have 99.7 percent of the theoretical copper content.

\*\*These values of  $T_s$  are derived from the values of  $A_s$ .

$$A_s = -\log_{10} T_s$$

Table 5. Values of Spectral Absorbancy,  $A_s$ , and Transmittancy,  $T_s$  for Standard Cobalt Ammonium Sulphate Solution, Unfiltered, Thickness 10.00 mm, Temperature 25°C having the following composition:\*

Cobalt ammonium sulphate ( $\text{CoSO}_4 \cdot (\text{NH}_4)_2\text{SO}_4 \cdot 6\text{H}_2\text{O}$ )\*\* 14.481 grams  
Sulphuric acid (specific gravity 1.835) 10.0 cc  
Water (distilled) to make 1000. cc

Wavelength mμ	$A_s$	$T_s$ ***	Wavelength mμ	$A_s$	$T_s$ ***
350	0.0038	0.991	600	0.0137	0.969
60	.0040	.991	10	.0124	.972
70	.0050	.989	20	.0115	.974
80	.0065	.985	30	.0112	.975
90	.0088	.980	40	.0110	.975
400	.0125	.972	650	.0105	.976
10	.0168	.962	60	.0097	.978
20	.0224	.950	70	.0087	.980
30	.0340	.925	80	.0076	.983
40	.0522	.887	90	.0066	.985
450	.0773	.837	700	.0054	.988
60	.1031	.789	10	.0046	.989
70	.1213	.756	20	.0038	.991
80	.1349	.733	30	.0032	.993
90	.1472	.713	40	.0030	.993
500	.1635	.686	750	.0028	.994
10	.1742	.670			
20	.1689	.678	Hg 404.7	.0144	.967
30	.1452	.716	Hg 435.8	.0437	.904
40	.1113	.774	Hg 491.6	.1497	.708
			He 501.6	.1661	.682
550	.0775	.837			
60	.0496	.892	Hg 546.1	.0901	.813
70	.0308	.932	Hg 578.0	.0219	.951
80	.0207	.953	He 587.6	.0167	.962
90	.0158	.964	He 667.8	.0089	.980

\* These data apply accurately also from 400 to 750 mμ to a similar solution made up with 10.3 grams of cobalt sulphate ( $\text{CoSO}_4 \cdot 7\text{H}_2\text{O}$ ).

\*\*Chemical analysis showed a cobalt (plus nickel) content of 100.0 percent of the theoretical, the ratio of nickel to cobalt (metals) being 1 to 200.

\*\*\*These values of  $T_s$  are derived from the values of  $A_s$ .  
 $A_s = -\log_{10} T_s$ .

The published values for absorbancy and transmittancy are given in tables 4 and 5 herein, together with the composition and certain other pertinent information. Many additional details are given in M114 (15) including the changes in absorbancy ( $A_s$ ) with temperature. Both solutions obey Beer's law over a considerable range of concentrations. The values given for the eight Hg and He wavelengths are considered the most reliable, with an uncertainty in  $A_s$  not exceeding 0.001 for the particular chemicals used. Spectrophotometric reproducibility of the chemicals is also considered in the paper. By increasing the thickness or concentration a wide range of the photometric scale can be covered, except at the shorter wavelengths.

2. An aqueous solution of potassium chromate,  $K_2CrO_4$  (0.04g/l) in 0.05N KOH. This solution has been studied by many investigators here and abroad, and is considered one of the most suitable as a standard of spectral transmittancy and absorbancy in the ultraviolet. The absorption in the violet is also useful because the copper and cobalt solutions have too little absorption in this region to be of much value. The most recent determination is given in NBS Research Paper 2331 (16), which also notes most of the previous work. This solution has also been used in a recent comparative survey of photoelectric spectrophotometers (17). The values adopted in RP2331 are given in table 6. These data were obtained from solutions prepared in the following ways:

(1) A solution of  $K_2CrO_4$ , stock material, reagent grade, 0.0400 gram per liter, in 0.05N KOH.

(2) A solution of  $K_2CrO_4$  of the same concentration and alkalinity as (1) but prepared from 0.0303 gram of  $K_2Cr_2O_7$ , which when converted gave 0.0400 gram of  $K_2CrO_4$  per liter.

At wavelengths greater than 260  $m\mu$ , both solutions were found to remain stable (in transmittancy) for 5 or 6 years when stored in ordinary storeroom glass bottles. "Flaking" may occur during this time and any sediment should be allowed to settle to the bottom of the bottle. It is recommended that alkali-resistant ware, now available, be used for storing the solutions. At wavelengths less than 260  $m\mu$ , it is recommended that solutions not over 6 months old prepared from chemicals of the highest purity be used.

Table 6. Values of Spectral Transmittancy,  $T_s$ , and Absorbancy,  $A_s$ , for Standard Potassium Chromate Solution, Unfiltered, Thickness 10.00 mm, Temperature 25°C, having the following composition:

0.0400 gram per liter of potassium chromate ( $K_2CrO_4$ ) in 0.05 normal potassium hydroxide solution\*

Wavelength m $\mu$	$T_s$	$A_s^{**}$	Wavelength m $\mu$	$T_s$	$A_s^{**}$
220	0.358	0.446	335	0.600	0.222
25	.601	.221	40	.483	.316
30	.674	.171	45	.373	.428
			50	.276	.559
35	.616	.210	55	.199	.701
40	.507	.295			
45	.402	.396	60	.148	.830
50	.319	.496	65	.116	.936
53.6	.279	.554	70	.103	.987
			75	.102	.991
55	.268	.572	80	.117	.932
60	.233	.633			
65	.202	.695	85	.150	.824
70	.180	.745	90	.202	.695
75	.175	.757	95	.294	.532
			400	.402	.396
80	.194	.712	04.7	.515	.288
85	.257	.590			
90	.373	.428	10	.632	.199
95	.533	.273	20	.751	.124
96.7	.598	.223	30	.824	.084
			35.8	.859	.066
300	.709	.149	40	.882	.054
02.2	.771	.113			
05	.834	.079	50	.927	.033
10	.895	.048	60	.960	.018
13.2	.905	.043	70	.980	.009
			80	.991	.004
15	.900	.046	90	.997	.001
20	.864	.064	500	1.000	.000
25	.804	.095			
30	.710	.149			
34.2	.620	.208			

\* This solution of potassium hydroxide can be prepared with sufficient accuracy by dissolving 3.3 grams of potassium hydroxide sticks (85% KOH) of reagent quality in sufficient distilled water to make 1 liter.

\*\*These values of  $A_s$  are derived from the values of  $T_s$ ,  
 $A_s = -\log_{10} T_s$

Note: Distilled water only was used in the reference cell.

The data of table 6 are based on extensive measurements made with the Hilger sector-photometer photographic spectrophotometer, and the Beckman DU photoelectric spectrophotometer, supplemented with data obtained with the Brackett quartz photoelectric spectrophotometer and, above 400  $m\mu$ , with data obtained on the König-Martens visual spectrophotometer and the General Electric recording spectrophotometer. The temperatures were kept close to 25°C. In the values of  $T_s$  given in table 6 there is considerable uncertainty in the third decimal.

#### 4. Reference Standards of Spectral Directional Reflectance Calibrated Relative to Freshly Prepared Magnesium Oxide.

As a fundamental standard of spectral directional reflectance nothing has as yet been found more suitable than freshly prepared magnesium oxide. Its (total) luminous reflectance is high, 0.97 or 0.98, and nothing has been found of certainly higher reflectance. Its luminous directional reflectance,  $R_{0.45}$ , is 1.00, and its spectral selectivity throughout the visible spectrum appears to be less than 1 percent. These data are based on work by Priest (18), McNicholas (19), and Preston (20) and are summarized in National Bureau of Standards letter circular LC-547 (21). More recent work by Benford and others (22, 23) and by Middleton and Sanders (24, 25) closely confirm these results and extend the data into the ultraviolet and infrared.

Other agencies besides the Bureau have also recommended MgO for the same or similar purpose. In 1931, the International Commission on Illumination adopted a resolution which may be translated as follows: "For the colorimetric measurement of opaque materials the luminance of the specimen studied ought to be expressed as a function of the luminance of a surface of the oxide of magnesium considered under the same conditions of illumination and observation" (26). In 1944, this method was incorporated in ASTM Standard Method of test for spectral characteristics and color of objects and materials (27).

While the characteristics above noted make fresh MgO an excellent fundamental standard, it has other characteristics that are undesirable and that make the calibration and use of a secondary working standard a very advisable procedure. An MgO surface is extremely fragile and thus

is not very suitable for continued handling. A more serious defect is that its spectral reflectance may change by slight but definite amounts within a short time (sometimes in a day) after preparation, the reflectance decreasing below 550 m $\mu$ . Furthermore the nature and extent of the changes seem somewhat variable (1, Fig.33;24). A third reason for use of a working standard is that slight variations in reflectance (0.1 or 0.2 percent) have been noted for different preparations of freshly prepared MgO. By calibrating the working standard against several different MgO preparations a more representative standard is obtained than would be any single MgO surface by itself.

White structural glass by the name of Vitrolite, with one surface polished, has proved the most suitable for working standards of spectral directional reflectance, although the material is not uniformly good for this purpose and must be selected with care. A considerable supply of suitable Vitrolite has been obtained by the National Bureau of Standards, and standards are calibrated and issued for either the General Electric recording spectrophotometer or the Beckman DU spectrophotometer.

#### 4.1 For the General Electric Recording Spectrophotometer.

The Vitrolite standards issued for use with the General Electric recording spectrophotometer are about 4 inches square and 5/16 inch thick. They are covered in NBS test fee schedule 202.185, items; to K. ~~Item i (now obsolete) refers to data taken with an early model of the recording spectrophotometer where the radiant energy is incident in a slightly diverging beam whose axis is perpendicular to the sample. The specular component of the reflected energy from the polished surface thus returns towards the entrance aperture, part of it passing out through this aperture and being lost, and part being intercepted by the sphere and contributing to the measurements. The values obtained and reported are strictly valid on another spectrophotometer only if the same fraction of the specular component of the reflected energy is intercepted by the sphere.~~

In a later model of the spectrophotometer the design is such that the radiant energy is incident in a slightly diverging beam whose axis is at 6° to the perpendicular to the sample. The specular component of the reflected

energy is thus diverted away from the entrance aperture towards a port on the side. This port may be filled with MgO or with a black material, so that for plane glossy surfaces the measurement can be made with the specular component "included" or "excluded". This is covered in test fee items j and K.

To give one a better idea of the spectral directional reflectance relative to MgO for the Vitrolite standards thus issued, there is given in table 7 a set of values from 400 to 1080  $m\mu$  with specular component excluded that apply to one of the NBS standards. Values reported under test fee items j and K with specular component excluded will probably be closely similar to these. With specular component included, the values are greater by roughly 0.04. ~~Under item i, with specular component partly included and partly excluded, the values are of course intermediate between the other two sets.~~

Only one Vitrolite working standard is needed for the measurement of spectral directional reflectance on the General Electric recording spectrophotometer. This calibrated Vitrolite standard and the samples to be tested are in turn placed at the sample aperture of the integrating sphere, and any highly reflecting substance such a MgO or  $MgCO_3$  may be used at the comparison aperture provided the material to be tested does not reflect more than the comparison material.

The directional reflectances of the test samples relative to freshly prepared MgO are then obtained by use of correction factors derived from the ratios between the standard Vitrolite values and the values for the Vitrolite read from the curve sheet. This is explained in detail in the report accompanying each standard.

By this procedure only one Vitrolite standard is necessary and the Bureau does not issue these standards in pairs, as some have requested. For transmittance measurements any two nearly identical white surfaces are suitable and no standard reflecting surface is necessary.

Table 7. Spectral Directional Reflectance of NBS Standard Vitrolite V1-G4 Relative to Freshly Prepared Magnesium Oxide for Excluded Specular Component of Reflected Radiant Energy on a General Electric Recording Spectrophotometer.

For Visible Spectrum (400 to 750 $m\mu$ ) (black velvet port) (10 $m\mu$ slits)		For Near Infrared Spectrum (730 to 1080 $m\mu$ ) (black cavity port) (20 $m\mu$ slits)	
Wavelength ( $m\mu$ )		Wavelength ( $m\mu$ )	
400	0.875	730	0.861
10	.865	40	.858
20	.859		
30	.858	750	.856
40	.856	60	.854
		70	.851
450	.863	80	.849
60	.871	90	.848
70	.874		
80	.876	800	.845
90	.878	10	.843
		20	.841
500	.880	30	.839
10	.883	40	.836
20	.885		
30	.887	850	.834
40	.888	60	.832
		70	.830
550	.888	80	.828
60	.888	90	.826
	.887		
80	.885	900	.824
90	.883	10	.822
		20	.821
600	.880	30	.819
10	.877	40	.817
20	.874		
30	.871	950	.817
40	.869	60	.816
		70	.814
650	.868	80	.814
60	.867	90	.814
70	.867		
80	.866	1000	.813
90	.866	10	.812
		20	.812
700	.865	30	.811
10	.865	40	.811
20	.864		
30	.863	1050	.810
40	.862	60	.810
		70	.810
750	.862	80	.810

4.2 For the Beckman Quartz (DU, Non-Recording)  
Spectrophotometer.

The Vitrolite standards issued for use with the Beckman DU spectrophotometer are about 1 1/2 by 2 inches and 5/16 inch thick. They are covered in NBS test fee schedule 202.185, items 1, m and n. In this case the radiant energy is incident upon the sample in a nearly parallel beam whose axis is perpendicular to the sample. The reflected energy accepted for measurement is taken in an annular "beam" whose axis is closely 45° to the perpendicular in all directions, but the parts of which may vary in direction roughly from 35° to 55°. The specular component is thus excluded from the measurements and the values obtained and reported do not differ greatly from those shown in table 7.\* At 400 mμ and below, these values as obtained on NBS standard V2-B4 for the Beckman spectrophotometer are:

350 mμ	0.753
360	.812
370	.826
380	.809
390	.852
400	.868

It is important to note that the values of spectral directional reflectance obtained and reported for use on the Beckman DU spectrophotometer are not valid for use on the General Electric spectrophotometer, and vice versa, because of the notably different irradiation-reception conditions of the two types of instrument. Similarly the values reported for either instrument should not be used for other types of instrument, unless the geometrical conditions are sufficiently similar as to make the values valid for such purpose.\*

---

\*

In any comparison of results of directional reflectance measurements on instruments with differing geometries of irradiation and reception, it should be remembered that the directional reflectance of MgO also varies with the geometrical conditions. See first paragraph of Section 4, above.

---

## 5. References.

1. Kasson S. Gibson, Spectrophotometry (200 to 1000 millimicrons) NBS Circular 484 (1949).
2. Massachusetts Institute of Technology Wavelength Tables, measured and compiled under the direction of George R. Harrison (John Wiley & Sons, Inc., New York, 1939).
3. Spectrophotometry, report of O.S.A. Progress Committee for 1922-23, K. S. Gibson, Chairman, J. Opt. Soc. Am. & Rev. Sci. Inst. 10, 169 (1925).
4. N. C. Beese, Cesium vapor lamps, J. Opt. Soc. Am. 36, 555 (1946).
5. Kasson S. Gibson and Margaret M. Balcom, Transmission measurements with the Beckman quartz spectrophotometer, J. Research NBS 38, 601 (1947) RP1798; also J. Opt. Soc. Am. 37, 593 (1947).
6. M. Garbuny, Relative Intensitäten und Übergangswahrscheinlichkeiten der rotgelben Neonlinien, Zeit. Physik 107, 362 (1937).
7. H. J. McNicholas, Equipment for routine spectral transmission and reflection measurements, BS J. Research 1, 793 (1928) RP30.
8. K. S. Gibson, Direct-reading photoelectric measurement of spectral transmission, J. Opt. Soc. Am. and Rev. Sci. Inst. 7, 693 (1923); also described in J. Opt. Soc. Am. 21, 564 (1931).
9. H. H. Cary and Arnold O. Beckman, A quartz photoelectric spectrophotometer, J. Opt. Soc. Am. 31, 682 (1941).
10. Kasson S. Gibson and Harry J. Keegan, Calibration and operation of the General Electric recording spectrophotometer of the National Bureau of Standards, J. Opt. Soc. Am. 28, 372 (1938).
11. George L. Buc and E. I. Stearns, Uses of retardation plates in spectrophotometry, II, Calibration of wavelength, J. Opt. Soc. Am. 35, 465 (1945).

12. Harry J. Keegan and Kasson S. Gibson, On the use of working standards of didymium and Vitrolite glasses for spectrophotometric measurements, J. Opt. Soc. Am. 34, 770 (1944).
13. K. S. Gibson, Geraldine K. Walker, and Mabel E. Brown, Filters for testing the reliability of spectrophotometers, J. Opt. Soc. Am. 24, 58 (1934).
14. Kasson S. Gibson and Marion A. Belknap, Permanence of glass standards of spectral transmittance, J. Research NBS 44, 463 (1950) RP2093; also (condensed) J. Opt. Soc. Am. 40, 435 (1950).
15. Raymond Davis and K. S. Gibson, Filters for the reproduction of sunlight and daylight and the determination of color temperature, B. S. Miscellaneous Publication No. 114 (1931) M114.
16. Geraldine W. Haupt, An alkaline solution of potassium chromate as a transmittancy standard in the ultraviolet, J. Research NBS 48, 414 (1952) RP2331; also (condensed) J. Opt. Soc. Am. 42, 441 (1952).
17. Wallace R. Brode, John H. Gould, James E. Whitney, and George M. Wyman, A comparative survey of spectrophotometers in the 210-760  $m\mu$  region, J. Opt. Soc. Am. 43, 862 (1953).
18. Irwin G. Priest and J. O. Riley, The selective reflectance of magnesium oxide, J. Opt. Soc. Am. 20, 156 (1930).
19. H. J. McNicholas, Absolute methods in reflectometry, BS J. Research 1, 29 (1928) RP3.
20. J. S. Preston, The reflection factor of magnesium oxide, Trans. Opt. Soc. (London) 31, 15 (1929-30).
21. Preparation and colorimetric properties of a magnesium-oxide reflectance standard, NBS letter circular LC-547 (1939).

22. Frank Benford, Gwen P. Lloyd, and Sally Schwarz, Coefficients of reflection of magnesium oxide and magnesium carbonate, J. Opt. Soc. Am. 38, 445 (1948).
23. Frank Benford, Sally Schwarz, and Gwen P. Lloyd, Coefficients of reflection in the ultra-violet of magnesium carbonate and oxide, J. Opt. Soc. Am. 38, 964 (1948).
24. W. E. K. Middleton and C. L. Sanders, The absolute spectral diffuse reflectance of magnesium oxide, J. Opt. Soc. Am., 41, 419 (1951).
25. C. L. Sanders and W. E. K. Middleton, The absolute spectral diffuse reflectance of magnesium oxide in the near infrared, J. Opt. Soc. Am. 43, 58 (1953).
26. Proceedings, C.I.E., 8th Session, Cambridge, p. 23, 1931.
27. A.S.T.M. standard method of test for spectral characteristics and color of objects and materials, ASTM D-307-44.

Letter Circular 1017  
(Supersedes LC 929 dated 11/26/48)  
Prepared by: KSG:HJK/EMR  
Photometry and Colorimetry Section  
Optics and Metrology Division  
National Bureau of Standards

Rerun October 1, 1962.